




Applanation Tonometer AT-1

Aiming at new levels in quality



Applanation Tonometer

AT-1

Applicable tonometer base

For SM-70N



S12-13

For SM-30N and SM-10N



S12-14

For SM-4N and SM-2N



S12-15

For 300XL




S12-16

Accessories

(Standard for all models)



Mounted on 300XL



Specifications

| Model | AT-1 | | | | |
|--|---|-----------------|------------------------------|--------------------------|--------------|
| Model of slit lamp | SM-90 SM-90N | SM-70 SM-70N | SM-30,SM-10 SM-30N,SM-10N | SM-4,SM-2 SM-4N,SM-2N | 300XL |
| Tonometer main unit | S12-12 | | | | S12-12 XL |
| Tonometer base | --- | S12-13 | S12-14 | S12-15 | S12-16 |
| Type of AT-1 | AT-1 | AT-1 Type 7 | AT-1 Type 3 | AT-1 Type 4 | AT-1 Type XL |
| Measurement method | Depression type (observed with left eye on the slit lamp) | | | | |
| Depression area | 7.354mm ² | | | | |
| Diameter of measurement depression area | 3.06mm | | | | |
| Measurement range | 0 ~ 80mmHg | | | | |
| Minimal scale | 2mmHg | | | | |
| Measurement error | ±1mmHg | | | | |
| Equation for conversion from mmHg to kPa | 1mmHg = 101325/760Pa ÷ 133.3Pa 100Pa = 0.1kPa | | | | |

● Design and specifications are subject to change as improvements are made to the product.

TAKAGI SEIKO CO.,LTD.

330-2 IWAFUNE, NAKANO-SHI, NAGANO-KEN, 383-8585, JAPAN

TEL.+81-269-22-4512 FAX.+81-269-26-6321

URL:<http://www.takagi-j.com> E-mail:info@takagi-j.com

