



KANGHUA



OPHTHALMIC EQUIPMENT PORTABLE SLIT LAMP MICROSCOPE

WE MAKE VISION MORE WONDERFUL

Using

1. Portable slit lamp with high-brightness, cold-light illuminator, without external power supply.
2. High-brightness white LED, cold-light illuminator, its brightness up to 30000 LUX, service life more than 50000 hours.
3. Sharp slit, the width can be adjusted continuously.
4. Telescopic forehead-rest, easier to observe and adjust the diopter.
5. Easier to adjust the magnification by removable objective lens.
6. Small and light, only 700g including battery.
7. Different filters: Heat Absorption; Yellow; Redfree (Green); Cobalt Blue.
8. Full Power supplied by battery can drive the slit lamp for 5-6 hours.



10X eyepiece



Diopter adjusting shaft



Forehead-rest



Box



Technical Specification

Item \ Model		SLM-6L	SLM-6M	SLM-6H
Microscope	Microscope Type	Converging		
	Objective Lens	1.0x	1.0x 1.6x	1.0x 1.6x
	Eyepiece	10x	10x	10x 16x
	Total Magnifications	10x	10x 16x	10x 16x 25.6x
	Working Distance	81mm		
	Intraocular Distance	45mm-78mm		
	Diopter Adjustment	-8D ~ +8D		
Illumination	Slit Width	0-12mm Continuous		
	Slit Height	1~12mm	1~12mm	1-12mm
	Slit Angle	Horizontal±30°		
	Light Spot Diameter	Φ 12mm	Φ0.2mm; Φ1mm; Φ5mm; Φ12mm	
	Filters	Heat Absorption; Yellow; Redfree (Green); Cobalt Blue		
	Illumination Bulb	High- Brightness White LED		
	Brightness	>30000Lux		
	Power	5VA		
Power	Charging Voltage	AC 100 ~240V, 50/60Hz		
	LED Working Voltage	DC 3V, 350mAh		
	Battery	Li-ion Battery; 7.4V, 680mAh. Full Battery can supply 5-6 hours continuously		
Weight and Dimension	Dimension	169.5 x 109.5 x 228.5mm		
	Weight	700g (Including Battery)		