

YURATEK

DURU

A professional Approach to
Refraction Units



www.yuratek.com

TECHNICAL SPECIFICATIONS UNIT AND CHAIR

Unit	DURU
Power	100-240 VAc 50-60 Hz
Electrical Classification	Class I, Type B
Table Top Electrical Output	110 or 220 VAc & 0-6/12 (S.L)
Head Lamp Dimmable	Led Lamp
Maximum Load of Table	20 Kg + 20 Kg + 20 Kg
Chair Stroke	180 mm

STANDARD FEATURES

Table top for 2 Instruments or 3 Instruments
Easy changes to Left or Right Hand Versions at site.
Electromagnetic Brake in Every Angle
Dimmable of Slit Lamp Control
Standard Trail Lens Drawer with self
Backlight Illuminated Keypad
Powered Elevation of Table Top (0-6/12V & 230 V Ac)
Knee Safety Switch
Energy Saving Mode
Mini Keypad on table

CONTROL PANEL

Power On/Off
Chair Up-Down Elevation Control
Aux On/Off Control
Projector On/Off Control
Head Lamp Dimmable Control
Slit Lamp Illumination Light Control
Table Up-Down Elevation Control
Motorized Phoropter Arm Control
Electro Magnetic Brake Control
Reading Lamp (Nearvision) Control

OPTIONAL ACCESSORIES

Wheel Chair Access Versions
Elevation Table Top
Inclinable Chair Backrest 45° or Full
Foot Switch for Electro-Brake
Manual Horizontal Chair Movement
LCD Monitor Support
Ophthalmoscope Retinascope Holder
Manual Phoropter Arm
Motorized Phoropter Arm

DURU



IMPORTANT Yuratek Ltd. subject to change in design and /or specifications without advanced notice for detailed info please contact your dealer.

YURATEK

YURATEK LTD. ŞTİ.
İvedik Org. San. Böl. 1475 Sk. No: 4 Ostim - ANKARA /TÜRKİYE
Tel : +90.312 395 29 96 Fax : +90.312 395 29 97
e-mail: info@yuratek.com www.yuratek.com

