

AJL
vet

LIOCAN CANINE PRELOADED IOL

PRODUCT DESCRIPTION AND INDICATION

Ophthalmic medical device designed for implantation in canine crystalline extraction procedures.

The product is made of HEMA (UV), a biocompatible acrylic copolymer that offers excellent geometric and folding properties allowing the lens implantation through a small incision.

Indicated for implantation in the posterior chamber, replacing the canine crystalline extracted for the correction of aphakia.

Veterinary use only.

HOW SUPPLIED

Single and sterile, preloaded in a disposable injector, ready to use.



PRODUCT SPECIFICATIONS AND MODELS

Optic Design

Biconvex
6.5 mm diameter
+41 D

Available Lengths and References

10.25 mm	[ref. TXY651025PL]
11.25 mm	[ref. TXY651125PL]
12.25 mm	[ref. TXY651225PL]
13.25 mm	[ref. TXY651325PL]
14.25 mm	[ref. TXY651425PL]